

# ENDODRY™

Storage & Drying System | **Dry & Store**





## THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, only Cantel Medical delivers the Complete Circle of Protection, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimize your success.

## A Higher Standard for Endoscope Drying and Storage

### Ensure Your Flexible Endoscopes Remain Safe and Dry During Storage

MEDIVATORS Storage and Drying Cabinet continuously forces HEPA filtered air through all endoscope channels, ensuring the endoscope is dry and fluid-free.

- Continuous ventilation eliminates residual moisture that can lead to bacteria growth
- Circulation of dry HEPA filtered air through each endoscope channel and around the outer sheath ensures the endoscope is dry prior to use.
- Short drying cycle improves scope turn-around time

### Protect Flexible Endoscopes from Damage

Cassette system eliminates free swinging of the distal tip and protects the endoscope during transport and storage.

- Transporting flexible endoscopes in cassettes reduces endoscope damage and related repairs by 25-30%
- Quick, efficient transfer from automated endoscope reprocessor to storage cabinet via cassette system

### Customizable to Fit Your Facility's Needs

Space-conscious design available in various configurations to meet the requirements of any facility.

- Available in single-sided or pass-thru models (front and back door) for separation of wet and dry areas
- Compatible with ADVANTAGE PLUS™ Reprocessor or DSD hook-ups to increase department efficiencies

### Maximize Patient Safety

MEDIVATORS Storage and Drying Cabinet uses a cassette to minimize touching and potential recontamination.

- Reduced worker handling minimizes endoscope recontamination and optimizes infection control
- No disconnection and reconnection of endoscope channels to hook-up ensures proper connectivity and time savings
- SMART Light (green, blue, red) for quick cabinet status identification

### Data Management Ensures Endoscope Drying Cycle Compliance

Continuous data monitoring provides detailed information on the drying status of all endoscopes.

- Constant air circulation controls temperature and humidity within the cabinet
- Barcode scanner for ease of endoscope ID and endoscope tray location; touchscreen for quick visualization
- SMART Connect data interface with ADVANTAGE PLUS Reprocessor and ENDORA™ Endoscope Tracking System



## DRY & STORE

Bacteria pose significant risk to endoscopes during endoscope storage. Cantel's transport, drying and storage solutions are designed to protect valuable inventory, reduce cross contamination touchpoints, eliminate moisture in the endoscope channels and control humidity. Higher humidity or moisture in endoscope channels is known to aid bacterial growth.

# ENDODRY™

Storage & Drying System | **Dry & Store**

Color touch screen display for quick identification of endoscope storage time



SMART Connect to ADVANTAGE PLUS™ Reprocessor and ENDORA™ Tracking Systems

Locking doors protect endoscope inventory

Available in left or right door swing configurations



Barcode scanner for ease of endoscope and operator identification



Compatible with ADVANTAGE PLUS Reprocessor or DSD hook-ups to increase efficiencies in your department



Convenient transport cassette for endoscope protection and recontamination prevention

Pass thru and single sided models available



Wet Side

Dry Side

## ENDODRY™ Cabinet Specifications

### Electrical Requirements

100-240 VAC, 50/60HZ, Power Input: 75W,  
Fuse: T2A (internal device) 15A (circuit)

### Electrical Safety Certifications

IEC 61010-1:2010

### Operating Temperature

50 - 104°F (10 - 40°C)

### Interface

RJ45 (10/100 Mbit)

### Dimensions

75 H x 24 W x 22 D (inches)

190 H x 60 W x 54 D cm

Door Swings 20.5 (inches) 52 (cm)

### Weight

375 lbs (170 kg)

### Endoscope Traceability

ENDORA™ Tracking System Integration

### Air Requirements

(Refer to technical documents)

Instrument grade air

## ENDODRY™ Cabinet Ordering Information

Model Number	Description
DRY-2001, DRY 2002, DRY-2003, DRY-2004**	ADVANTAGE PLUS™ AER – Pass-thru model
DRY-1001, DRY-1002*	ADVANTAGE PLUS AER – Single-side entry
DRY-2005, DRY 2006, DRY-2007, DRY-2008*	DSD AER Platform – Pass-thru model
DRY-1003, DRY-1004*	DSD AER Platform – Single-side entry
CAS-1000	Endoscope Transport Cassette

\*\* Available in right or left door swing configurations

\*Hook-ups sold separately.

ADVANTAGE PLUS™, ENDORA™, ENDODRY™ and MEDIVATORS™ are trademarks of Medivators Inc.



**Cantel Medical (Canada)**  
90 Gough Road, Unit 1  
Markham, ON L3R 5V5

Phone: +1.844.348.5636  
Fax: +1.844.348.5637  
Email: [info@cantelmedical.ca](mailto:info@cantelmedical.ca)

[www.cantelcanada.com](http://www.cantelcanada.com)

© 2017 Medivators Inc. P/N 50098-710D-ENCAN