

SCOPE BUDDY™

Endoscope Flushing Aid | **Manual Clean***



*SCOPE BUDDY is only intended to be used as a Flushing Aid during Manual Cleaning step of Endoscope Reprocessing. It is not intended to manually clean endoscopes.



THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, only Cantel Medical delivers the Complete Circle of Protection, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimize your success.

Consistent, Hands-free Flushing of Endoscopes

Reduces the Risk of Repetitive Motion Injuries

- Replaces repetitive and labor-intensive syringing of channels
- Simple, hands-free operation reduces muscle fatigue
- Provides constant fluid delivery through all endoscope channels

Provides Consistent Flushing

- Continuous, timed fluid delivery ensures consistent, repeatable endoscope flushing
- Audible and visual alerts upon completion
- Compatible with all major endoscope brands
- Works with your existing sink; compact size for countertop or wall mounting
- Meets or exceeds endoscope manufacturers requirements for channel flushing volumes

Shortens Flushing Time

- Thoroughly flushes most endoscopes in 35 seconds
- Flushes multiple endoscope channels simultaneously
- Quickly and quietly delivers fluids through endoscope channels

INTERCEPT™ Detergent Provides Superior Endoscope Cleaning

- Fast, one-minute contact time
- Low foaming, neutral pH formulation safe for endoscope use
- Cost-effective concentrated formula reduces per-use costs



MANUAL CLEAN



Manual cleaning is essential in removing bioburden and preparing the endoscope for high-level disinfection. For optimal results, Cantel's manual cleaning products and detergents are designed to be a critical first step to efficiently and effectively eliminating bacteria on endoscopes prior to reprocessing.

SCOPE BUDDY™

Endoscope Flushing Aid | **Manual Clean***



Visual display of endoscope flushing time

User-selectable time setting

Audible and visual alerts upon completion

One-touch start/stop operation

Simple and quick tubing connections

*SCOPE BUDDY is only intended to be used as a Flushing Aid during Manual Cleaning step of Endoscope Reprocessing. It is not intended to manually clean endoscopes.

Complete Endoscope Connection Kit



Wall Bracket

Power Supply

SCOPE BUDDY™ Flushing Aid

Unit Dimensions

9.5 H x 6 W x 8 D (inches)
 (24.1 H x 15.2 W x 20.3 D (cm))

Unit Weight

10 lbs (4.5 kg)

Electrical Requirements

110 – 230 V, 50/60 Hz

Electrical Safety Certifications

UL Standard 61010-1
 CANADA/USA C22.2 No. 61010-1
 ETL and cETL listed

Fluid Delivery Rate

1 liter/minute

Flush Time

0 – 99 minutes

Timer Control

Digital display with touchpad controls

System Pressure

Built-in safety bypass

Cycle Completion Indicators

Audible and visual alerts

Unit Placement

Counter-top or wall-mounting
 (bracket included)

Warranty

One year



SCOPE BUDDY Flushing Aid Ordering Information

Model Number	Description
EFA-US-G	North America
EFA-IN-G	International

Accessory Ordering Information

Part Number	Description
78399-669	Complete Endoscope Connection Kit
75227-679	Graduated Cylinder

Detergent Ordering Information

Part Number	Description
ML02-0133	INTERCEPT™ Detergent (4 one-gallon bottles/case)

SCOPE BUDDY™ and INTERCEPT™ are trademarks of Medivators Inc.



Cantel Medical (Canada)
 90 Gough Road, Unit 1
 Markham, ON L3R 5V5

Phone: +1.844.348.5636
 Fax: +1.844.348.5637
 Email: info@cantelmedical.ca

www.cantelcanada.com

© 2017 Medivators Inc. P/N 50098-077C-ENCAN