



Secure. Tested. Proven.

We're your partner for sterilization products.



Monitoring Products for Steam Sterilization Processes

STEAM

Select a system that works in all steam sterilization cycles. Cantel Medical products for monitoring steam sterilization processes provide fast, visible, consistent and reliable results offering peace of mind for assurance of effective sterilization.

STEAMPlus™ Sterilization Integrators

- Immediate read-out with wicking-style movement
- Reacts to all three parameters of sterilization (time, temperature, saturated steam)
- Larger size four inch strip has distinct pass/fail criteria
- ISO 11140-1 Type 5 Compliant
- Lot number and expiration date displayed on front for easy reading



Bowie-Dick Type Test Packs

- Designed to detect the presence of air in prevacuum steam sterilizers
- Available in single-use or reusable test packs
- Easy-to-interpret indicator sheet with blue-to-dark color transition
- Instructions for use printed on the back of each pack
- ISO 11140-5 Type 2 Compliant

SporView® Self-Contained Biological Indicators

- Compatible with all steam sterilization systems
- *Geobacillus stearothermophilus* product configurations available
- Visually distinct positive results from purple to yellow
- Chemical Indicator on label provides immediate indication of exposure
- Spore growth of a test BI is observable in as little as 2 ½ hours with 10-hour final incubation
- Compliant with ISO 11138-1 and ISO 11138-3
- Available as a single use test pack with Type 5 integrator for immediate readout
- Media growth promotion has been validated for growth of *G. stearothermophilus* organism after exposure to extended steam cycles



Steam Sterilization Pouches

- Single-use pouches come in self-seal, heat-seal and tubing for custom length
- Manufactured with medical grade paper, extra wide chevron seal, tacked corner seals and thumb notch for ease of opening
- Sure-Check pouches include internal/external multi-parameter indicators that respond to all 3 sterilization variables
- Available in multiple sizes
- Meets ISO 11607 standards for pouch performance
- Validated for double pouching*



Universal Lot Control Labels

- Manufactured free of natural rubber latex
- Multiple styles and colors available
- Can be used in steam, EO gas and hydrogen peroxide cycles

Indicator Tapes

- Designed to aid in distinguishing processed from unprocessed packs
- Consists of a strong custom-made white polyester substrate with a pressure sensitive adhesive which adheres to plastic, woven, treated woven and non-woven wraps
- Easy to remove adhesive allows for secure seals
- LF2 tape is manufactured free of latex & lead
- ISO 11140-1 Type 1 Compliant



Chemical Indicator Strips

- Certified to meet chemical indicator specifications for pack monitoring when used in steam sterilization processes
- Provides a distinct color change upon exposure to steam processes
- Demonstrates exceptional performance in the most widely used steam sterilizers
- ISO 11140-1 Type 4 Compliant

PRODUCT IMAGE	DESCRIPTION	CATALOG NO.†	UOM	SOLD PER
	STEAMPlus™ Sterilization Integrator (Type 5)	SSI-100	100/Pack	Pack
		SSI-1000 (4")	1000/Pack	Pack
	SSI-1000E (8")			
	STEAMPlus™ Test Pack	STP-025	24/Case	Case
	Tray Record Card	TRC-050	50/Box	Box
	Steam Indicator Strip (Type 4)	SIS-250 (4")	250/Box	Box
		SIL-250 (8")		
	Steam/EO Indicator Strip (Type 1)	SIL-1000 (Circle)	1000/Roll, 10 Roll/Case	Roll
	AirView™ Bowie Dick Test Pack	MBD030	30/Case	Case
		SBD-030		
	Reusable Lantor Cube® Test Pack	LCR-025	25/Cube	Cube
	Reusable Lantor Cube® Metal Clamp	LCP-001	1 Ea.	Ea.
	SporView® Self-Contained Biological Indicator	SVT-050 (10 hours)	50 Vials/Box	Box
		SCS-025 (24 hours)	25 Vials/Box	
		SCS-100 (24 hours)	100 Vials/Box	
	SporView® Plus Test Pack	SBT-025, SBT-255	25/Case	Case
	Sure-Check® & Duo-Check® Sterilization Pouches	SC-Series 2"-18"	100-200/Box, 5-20 Boxes/Case	Case
	Sani-Roll® Sterilization Tubing	SCT-Series 2"-10"	100ft/Roll, 6-16 Rolls/Case	Case
	Universal Lot Control Label	URL-series	1125/Roll, 12 Rolls/Box	Box
	Chemical Indicator Tape	ST-036, ST-048, ST-072	1/Roll	Roll
	Blue Steam Indicator Tape	BT-036, BT-048		
	Lead/Latex Free Indicator Tape	LF2-036, LF2-048		

Protect with the pros you can trust: Cantel

A leading global company dedicated to delivering innovative infection prevention and control and sterility assurance solutions, Cantel offers a broad line of high-quality products for monitoring steam sterilization. Each is backed by a team of dedicated experts you can count on to be there to help you with compliance – how, when and where you need them – as you work to improve outcomes, enhance safety – and save lives.

* Study on file † Additional products available



Unless otherwise specified, all trademarks are property of Crosstex International, Inc or SPS Medical Supply Corp. Trypticase® is a registered trademark of Becton Dickinson and Company. Lantor Cube® is a registered trademark of Lantor UK LTD.

Cantel (Canada)
90 Gough Road, Unit 1
Markham, ON L3R 5V5

Toll Free Phone: +1.844.348.5637
Toll Free Fax: +1.844.348.5637
Email: info@cantelmedical.ca

www.cantelcanada.com