

# ENDOSCOPE

Care and Accessories | **Procedure**





## THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, only Cantel Medical delivers the Complete Circle of Protection, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimize your success.

## Simplify and Improve Specimen Retrieval During Colonoscopy

### Polyp Retrieval

Installing THE **ADVANTAGE** POLYP TRAP™ Polyp Capture System in-line between your endoscope and the suction chamber makes it fast and easy to recover patient polyps during your procedure. In-line capture facilitates uninterrupted suction and eliminates the need for personal handling, reducing the risk of lost tissue. Single-use containers safeguard your patients and increase infection prevention.

- Available in 1, 2 or 4-chambered models
- Features individually numbered and removable tissue chambers
- Packaged for single-patient use
- Compatible with all major endoscope brands
- Clear plastic construction for enhanced visualization potential cross contamination



ECONOTRAP™ Polyp Trap  
Single Chamber Polyp Trap



ECONOTRAP™ Polyp Trap  
Two Chamber Polyp Trap



THE **ADVANTAGE** POLYP TRAP™  
Collection Container (Four Chamber)



## PROCEDURE

Reducing the risk of patient cross contamination is at the forefront of infection prevention. Cantel innovates infection control products designed to improve patient outcomes, while increasing procedural efficiency.


# ENDOSCOPE

Care and Accessories | **Procedure**

## Adult and Pediatric Bite Blocks

Disposable Bite Blocks help protect both your patient and your equipment. Side ports allow you to reposition the endoscope, and facilitate access for suction.

Designed with the highest standards of patient comfort and safety in mind, our bite blocks offer a secure, comfort-fit design and an insertion port to accommodate a wide range of scope types. MEDIVATORS Bite blocks and straps are latex-free, remain in position during procedure and will not cause discomfort to your patients.



## Cleaning Brushes

MEDIVATORS cleaning brush line combines maximum cleaning capabilities with versatile brush options to meet your challenging reprocessing demands.

Both the single ended brush and the double ended brush offer desired catheter stiffness for ease of use and nylon bristles to provide the highest protection against channel damage.

### Polyp Trap Ordering Information

Part Number	Description	Units/Box
100101	ECONOTRAP™ 1 Polyp Trap (Single Chamber)	50 per box
100102	ECONOTRAP 2 Polyp Trap (Two Chamber)	50 per box
100104	THE <i>ADVANTAGE</i> POLYP TRAP™ Collection Container (Four Chamber)	50 per box

### Bite Block Ordering Information

Part Number	Description	Units/Box
100411	Adult Bite Block - 60 Fr Port	50 per box
100412	Pediatric Bite Block - 45 Fr Port	50 per box

### Cleaning Brushes Ordering Information

Part Number	Description	Units/Box
100401	Single Ended Cleaning Brush	50 per box
100402	Double Ended Cleaning Brush	50 per box

THE ADVANTAGE POLYP TRAP™ and ECONOTRAP™ are trademarks of Medivators Inc.



**Cantel Medical (Canada)**  
90 Gough Road, Unit 1  
Markham, ON L3R 5V5

Phone: +1.844.348.5636  
Fax: +1.844.348.5637  
Email: [info@cantelmedical.ca](mailto:info@cantelmedical.ca)

[www.cantelcanada.com](http://www.cantelcanada.com)

© 2017 Medivators Inc. P/N 50098-1407A-ENCAN