

VistaClear™

Centralized Water Filtration System

VistaClear.com

 ordervista@crostex.com  800-531-2221  vrg.support

Powerful protection and convenience, centralized.

The VistaClear™ Centralized Water Filtration System provides dental practices with centrally installed water filtration, backflow prevention and a method for convenient waterline cleaning — eliminating the need for dental bottles altogether! This patented, award-winning system saves practices time and money while protecting dental equipment.

CVV1000-28-HP



Saves time.

The VistaClear Centralized Water Filtration System's on-demand submicron water filtration eliminates wait time and the need for storage tanks and bottle systems. And its centralized, patented cleaning method makes periodic waterline maintenance easier and more efficient.



Saves money.

The VistaClear Centralized Water Filtration System includes CSA-certified backflow prevention and eliminates the need for reduced pressure zone (RPZ) devices and the ongoing cost of annual inspections. It also eliminates the need for purchased water and daily maintenance products.



Protects equipment.

With submicron filtration and no harmful chemicals, the VistaClear Centralized Water Filtration System protects expensive handpieces and turbines from the wear-and-tear typical water supplies can inflict. The VistaClear Centralized Water Filtration System filters remove virtually all contaminants that often reside in city water supplies.

Upgrade for peak performance.

Do you own an older VistaClear™ Centralized Water Filtration System? You can upgrade your original system to the current design for better performance and lower cost of ownership — for approximately the same price as replacing your filters. Contact your favorite dental dealer for more information.

CVR9280 VistaClear™ Centralized Water Filtration System HP Retrofit System



Serve up to 10 operatories.

The VistaClear Centralized Water Filtration System comes standardly equipped in an eight-operator configuration. It's possible to connect up to 10 operatories to one VistaClear Centralized Water Filtration System with an available accessory manifold.



CVDM125-10C

VistaClear™ Centralized Water Filtration System 10-Op Distribution Manifold



24"

16"

5 7/8"

Source Water Input Specifications

Operating Temperature Range
45-100°F (7.22-37.78°C)

Operating Pressure Range
30-100 psi (2.07-6.89 bar)

Optimal pH Range
6.5 – 8.5

Maximum pH Range
5.5 – 9.5

Maximum Influent Manganese
0.05 ppm

Maximum Influent Hydrogen Sulfide
0.5 ppm

Maximum Influent Chlorine / Chloramine
5 ppm

Maximum Continuous Flow Rate / Filter
0.2 gpm (0.76 Lpm)

Additional Information

Filter Element Service Life
12 months

Replacement Filter Kit
R9720

Dimensions (W x H x D)
16" x 24" x 5-7/8"

Approximate Ship Weight
28 lbs.

i Use of a VistaClear Centralized Water Filtration System does not, by itself, eliminate the need for proper user maintenance. Due to the nature and complex design of dental delivery units, the periodic use of waterline cleaners and antimicrobials is extremely important for proper dental waterline care. For this reason, Crosstex International, Inc. does not make specific bacteria reduction claims for the output water from dental delivery units connected to the VistaClear Centralized Water Filtration System.



All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.

© 2018 Crosstex International, Inc. 0119 DLIT00645 Rev B